
Question 4



Discussion:

The topomaps show a noticeable reduction in alpha power (8-12 Hz) from the baseline to the task condition for both subjects, particularly in the rear regions of the scalp. During the baseline, both Subject 06 and Subject 07 exhibit strong alpha activity, especially in the rear areas, which is typical for relaxed or resting states. In the task condition, both subjects show a decrease in alpha power, indicating a shift to a more focused and engaged mental state. This reduction, known as alpha desynchronization, is more pronounced and widespread in Subject 06 compared to Subject 07, suggesting greater cognitive engagement for Subject 06. Overall, the reduction in alpha power from baseline to task reflects the brain's transition from relaxation to active task performance.

Code:

```
import mne
```

```
import os
import numpy as np
import matplotlib.pyplot as plt
# Load subject 06 edf files
subject 06 baseline = mne.io.read raw edf(os.path.join("data", "Subject06 1.edf"),
preload=True)
subject 06 task = mne.io.read raw edf(os.path.join("data", "Subject06 2.edf"),
preload=True)
# Load subject 07 edf files
subject 07 baseline = mne.io.read raw edf(os.path.join("data", "Subject07 1.edf"),
preload=True)
subject 07 task = mne.io.read raw edf(os.path.join("data", "Subject07 2.edf"),
preload=True)
# Clean up channels
raw data = [subject 06 baseline, subject 06 task, subject 07 baseline,
subject 07 task]
for raw in raw data:
   raw.rename channels({
       'EEG Fp1': 'Fp1', 'EEG Fp2': 'Fp2', 'EEG F3': 'F3', 'EEG F4': 'F4',
       'EEG F7': 'F7', 'EEG F8': 'F8', 'EEG T3': 'T3', 'EEG T4': 'T4',
       'EEG C3': 'C3', 'EEG C4': 'C4', 'EEG T5': 'T5', 'EEG T6': 'T6',
       'EEG P3': 'P3', 'EEG P4': 'P4', 'EEG O1': 'O1', 'EEG O2': 'O2',
       'EEG Fz': 'Fz', 'EEG Cz': 'Cz', 'EEG Pz': 'Pz', 'EEG A2-A1': 'A2',
       'ECG ECG': 'ECG'
   })
   raw.set channel types({'ECG': 'ecg'})
  # Set the montage for the channel locations
  montage = mne.channels.make standard montage('standard 1020')
  raw.set_montage(montage)
alpha band = (8, 12)
# List of subjects and their corresponding data
subjects = [
   ('Subject 06 Baseline', subject_06_baseline),
```

```
('Subject 06 Task', subject_06_task),
   ('Subject 07 Baseline', subject 07 baseline),
   ('Subject 07 Task', subject 07 task)
fig, axes = plt.subplots(2, 2, figsize=(15, 12)) # 2x2 grid for 4 plots
# Loop through each subject and condition
for i, (title, raw) in enumerate(subjects):
  # Compute average PSD for each condition
  psd, freqs = raw.compute_psd(method='welch', fmin=alpha_band[0],
fmax=alpha band[1], n fft=2048).get data(return freqs=True)
   # Average the PSD across the frequency range for the alpha band
  alpha psd = np.mean(psd[:, (freqs >= alpha band[0]) & (freqs <= alpha band[1])],</pre>
axis=1)
  row, col = divmod(i, 2)
  # Create topomap
  im, cn = mne.viz.plot topomap(alpha psd, raw.info, axes=axes[row, col], show=False,
contours=8, cmap='viridis')
  axes[row, col].set_title(f'Alpha Band (8-12 Hz) PSD Topomap - {title}',
fontsize=14)
  axes[row, col].grid(True)
# Add color bar
cbar_ax = fig.add_axes([0.92, 0.15, 0.02, 0.7])
fig.colorbar(im, cax=cbar_ax, orientation='vertical', label='Power (dB)')
plt.tight layout(rect=[0, 0, 0.9, 1])
plt.show()
```